

LIFT HYGIENE CHAIR

Calypso™



An easy-to-use lift hygiene chair

Calypso is an integral part of a safe, comfortable bathing cycle. This lift hygiene chair from Arjo gives the semi-dependent resident constant security during transfers, transport and the bathing session itself.

With Calypso, a single carer can safely manage all aspects of the bathing cycle and work in an ergonomically sound way at all times. This battery-powered lift hygiene chair offers many features designed to enhance safe, ergonomic working practices.

Using Calypso in conjunction with a height-adjustable Arjo tub is the perfect combination for optimised ergonomics, as carers can always perform routines at the correct working height.

For more information visit www.arjo.com

arjo
with people in mind

Key features



BATHING

Calypso integrates with several height-adjustable Arjo bathing systems.



HANDLES

Easy-to-grip handles make manoeuvring and positioning smooth and easy.



TOILETING

With Calypso, manual lifting is eliminated from the routine and staff intrusion is kept to a minimum.



SEAT

Two seat versions are available: the Classic seat and a version in soft, easy-to-clean PUR material.



SAFETY BELT

The safety belt ensures resident security when seated and prevents the risk of sliding off the seat while bathing.



BEDPAN

The bedpan is inserted in the holder under the seat for the convenience of incontinent residents.

PRODUCT INFORMATION

Lifting capacity	136 kg (300 lbs.)
Weight	55 kg (120 lbs.)
Width	680 mm (26 3/4")
Length	815 mm (32")
Lifting stroke range	565 mm (22 1/4")
The seat is available in two versions:	Classic and in PUR material.
Powered functions	Raising & lowering
Bath compatibility	With Rhapsody/Primo (System 20-25), Rhapsody (System 25), Classic Line, Century
Best suited for Mobility levels	Albert, Barbara and Carl

All products are CE marked.

Only Arjo designed parts, which are designed specifically for the purpose, should be used on the equipment and products supplied by Arjo.
As our policy is one of continuous development we reserve the right to modify designs and specifications without prior notice.

® and ™ are trademarks belonging to the Arjo group of companies.
© Arjo, 2019

